

## Product datasheet for **RR201341L3V**

### **Snrpn (NM\_031117) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Snrpn (NM_031117) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Snrpn
Synonyms:	FE 294; FE294; FE 294 psi; FE294 psi
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031117
ORF Size:	720 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031117.2</a> , <a href="#">NP_112379.1</a>
RefSeq Size:	1428 bp
RefSeq ORF:	723 bp
Locus ID:	81781
UniProt ID:	<a href="#">P63164</a>
Cytogenetics:	1q36
Gene Summary:	component of the small nuclear ribonucleoprotein (snRNP) complex; may be involved in tissue-specific alternative RNA processing events, such as pre-mRNA splicing [RGD, Feb 2006]



[View online »](#)