

Product datasheet for **RR201355L3V**

Tfec (NM_022379) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tfec (NM_022379) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tfec
Synonyms:	Tcfec
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022379
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201355).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022379.1 , NP_071774.1
RefSeq Size:	1619 bp
RefSeq ORF:	954 bp
Locus ID:	26296
UniProt ID:	Q63302
Cytogenetics:	4q21
Gene Summary:	basic helix-loop-helix (BHLH) DNA-binding protein; contains a leucine zipper motif; forms heterodimers with TFE3 and inhibits TFE3-dependent transcription activation [RGD, Feb 2006]



[View online »](#)