

Product datasheet for **RR201383**

Arf3 (NM_080904) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arf3 (NM_080904) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Arf3
Synonyms: Ac1-253
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201383 representing NM_080904
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGCAATATCTTTGGGAACCTTCTGAAGAGCCTGATCGGGAAGAAGGAGATGCGCATCCTGATGGTGG
 GCCTGGATGCTGCAGGAAAGACCACCATCTCTATAAGCTGAAGCTTGGCGAGATTGTTACCACCATTCC
 CACCATCGGGTTAATGTGGAGACAGTGAATATAAGAATATCAGCTTACAGTATGGGACGTAGGTGGC
 CAGGACAAGATTCGGCCCTCTGGAGACTACTTCCAGAATACCAAGGCTTGATATTTGTGGTTGACA
 GCAATGACCGGGAGCGAGTGAACGAGGCCCGGGAAGAGCTGATGAGGATGCTGGCGGAGGATGAGCTCCG
 GGATGCTGTGCTGCTTGTGTTTGCAAACAAACAGGATTTGCCGAATGCTATGAATGCTGCTGAGATCACA
 GACAAGCTGGGCCTGCACTCTCTGCGCCACCGCAATTGGTACATTCAGGCCACCTGTGCCACCAGCGGGG
 ACGGGCTCTATGAAGGCTTGGACTGGCTGGCCAATCAGCTCAAAAACAAGAAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR201383 representing NM_080904
 Red=Cloning site Green=Tags(s)

MGNIFGNLLKSLIGKKEMRILMVGLDAAGKTTILYKLLKLGIVTTIPTIGFNVEVVEYKNISFTVWVDVGG
 QDKIRPLWRHYFQNTQGLIFVVDSDNRERVNEAREELMRMLAEDELRDVLLVFANKQDLPNAMAAEIT
 DKLGLHSLRHRNWIQATCATSGDGLYEGLDWLANQLKNKK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_080904

ORF Size: 543 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_080904.2</u> , <u>NP_543180.1</u>
RefSeq Size:	3315 bp
RefSeq ORF:	546 bp
Locus ID:	140940
UniProt ID:	<u>P61206</u>
Cytogenetics:	7q36
MW:	20.6 kDa
Gene Summary:	may play a role in intracellular vesicular transport [RGD, Feb 2006]