

Product datasheet for **RR201404L3V**

Bmp2 (NM_017178) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bmp2 (NM_017178) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bmp2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017178
ORF Size:	1179 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017178.1
RefSeq Size:	1275 bp
RefSeq ORF:	1182 bp
Locus ID:	29373
UniProt ID:	P49001
Cytogenetics:	3q36
Gene Summary:	involved in cellular signaling during limb development; induces bone formation [RGD, Feb 2006]



[View online »](#)