

Product datasheet for **RR201422**

Polr3k (NM_001014259) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Polr3k (NM_001014259) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Polr3k
Synonyms: MGC109430
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201422 representing NM_001014259
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTGCTGTTCTGCCCGGGCTGCGGGAACGGCCTCATTGTGGAGGAGGGGCAGCGCTGTCACCGCTTCG
CCTGCAACACGTGCCCGTACGTCCATAACATCACCCGCAAGGTAACAAACCGAAAGTATCCAAAGCTAAA
AGAAGTGGATGATGTGCTTGGTGGAGCAGCTGCCTGGGAGAATGTGGACTCTACTGCAGAGCCATGTCCA
AAATGTGAGCATCCTCGTGCCTACTTCATGCAGCTTCAGACCCGTTCTGCGGACGAGCCAATGACCACCT
TCTATAAGTGCTGCAATGCTCAGTGTGGACATCGCTGGAGGGATTCTATAAGTGCTGCAATGCTCAGTGT
GGACATCGCTGGAGGGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR201422 representing NM_001014259
Red=Cloning site Green=Tags(s)

MLLFCPGCGNGLIVEEGQRCHRFACTCPYVHNI TRKVTNRKYPKLKEVDDVLGGAAAWENV DSTAEPCP
KCEHPRAYFMQLQTRSADPEMPTTFYKCCNAQCGRWRDSISAAML SVDIAGG

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001014259

ORF Size: 324 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	995 bp
RefSeq ORF:	327 bp
Locus ID:	366277
Cytogenetics:	3q43
MW:	12.3 kDa
Gene Summary:	estrogen responsive protein; may be regulated by 17beta-estradiol [RGD, Feb 2006]