

Product datasheet for **RR201453L3V**

Capn10 (NM_031673) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Capn10 (NM_031673) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Capn10
Synonyms:	MGC93346
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031673
ORF Size:	1998 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201453).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031673.2
RefSeq Size:	2561 bp
RefSeq ORF:	2001 bp
Locus ID:	63834
UniProt ID:	Q9ES66
Cytogenetics:	9q36
Gene Summary:	large subunit of calpain, a calcium-dependent cysteine protease; human homolog has been associated with type II diabetes mellitus [RGD, Feb 2006]



[View online »](#)