

Product datasheet for **RR201472L4V**

Cxadr (NM_053570) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cxadr (NM_053570) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cxadr
Synonyms:	Car1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053570
ORF Size:	1056 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201472).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053570.1 , NP_446022.1
RefSeq Size:	1639 bp
RefSeq ORF:	1059 bp
Locus ID:	89843
UniProt ID:	Q9R066
Cytogenetics:	11q11
Gene Summary:	cell surface receptor for adenovirus and coxsackie virus; may be a barrier for macromolecule and ion movement across membranes in polarized epithelial cells [RGD, Feb 2006]



[View online »](#)