

Product datasheet for **RR201499L3V**

Pabpn1l (NM_001113779) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pabpn1l (NM_001113779) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pabpn1l
Synonyms:	Epabp2; Pabpn1; RGD1562761
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113779
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201499).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113779.1 , NP_001107251.1
RefSeq Size:	1154 bp
RefSeq ORF:	810 bp
Locus ID:	307920
UniProt ID:	B0BNE4
Cytogenetics:	19q12
Gene Summary:	Binds the poly(A) tail of mRNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)