

Product datasheet for **RR201517**

Ndufa11 (NM_212517) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ndufa11 (NM_212517) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ndufa11
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201517 representing NM_212517
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCAAGCGGTTTTTCGAAGCCTACAATGAGACCCCGACGGCACCCAGTGTATCGCAAGGCCTACA
TCACTACTGCCTTAGCGGCCTTTGTGGCCTAATAAGCTCTGCTTACAGCATCACACACAACCTGCAGA
TTCTCCCCTGGAAGCCGTGGCCAGGGTCGGCCGGTTCACATTCAGTGCAGCTGCTATTGGAGCTATGTTT
GGCCTTACCACCTGTGTGTCAGTGCCAGGTCCGAGAGAAGCCAGATGACCCCTGAACTACTTCATTGGAG
GTTGTGCGGGAGGCCTGACACTGGGAGCACGAACCTCACAGCTACGGGACCGCCGCAATTGCTGTGTCTA
CATGGGTACCGCGGCCGCCCTGTTCAAGATGGCAAGCTGGAGGGCTGGGAGCTGTTGCCACCCCAAG
GTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR201517 representing NM_212517
Red=Cloning site Green=Tags(s)
MAKRFFEAYNETPDGTQCHRKAYITTALGGLCGLISSAYSITHNPADSPLEAVARVGRFTFTAAAIIGAMF
GLTTCVSAQVREKPDPLNYFIGGCAGGLTLGARTHSYGTAAIGCVYMGTAALFKMGKLEGWELFATPK
V

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN:

NM_212517

ORF Size:

423 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_212517.1</u> , <u>NP_997682.1</u>
RefSeq Size:	538 bp
RefSeq ORF:	426 bp
Locus ID:	301123
UniProt ID:	<u>Q80W89</u>
Cytogenetics:	9q12
MW:	14.9 kDa
Gene Summary:	Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.[UniProtKB/Swiss-Prot Function]