

Product datasheet for **RR201521**

Nudt12 (NM_001109010) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nudt12 (NM_001109010) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nudt12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201521 representing NM_001109010
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTCTGTCAAAGAATCCCAAAGGAAATAATTTCTGAACTCCACAGTTCAGCCGAGAAGGAG
ACGTCGCGAAGTTAGCAGGAATTCAGTCATTCTCCATCTCTTCTCAACGAGACTTCTGAAAACGGATG
GACTGCTTTAATGTATGCCGCCAGGAACGGGCATCCGGATGCTGCCAGTTTCTGCTTGAGAAAGGATGT
GACAGATCGTCTGTCAATAAATCGAGGCAAACCGCCCTGGATATTGCGGTATTTGGGGTTATAAGCATA
TAGCTAATTTATTGGCTAACGCCAAAGGAGGGAAGAAGCCCTGGTTCTAACCAGTGACGTAGACGGATG
TGAAAATATTTTAGCAGAACTGCTGGACCGGAGGAGTGACAAAAGGACTAATTCTGACTGGCTGCAA
GCTAAAGAGAGCCACCCACCACAGTTTATATCCTTTTCTGACTTGAACCCCTGGTTACTCTAGGTG
GTAATAAAGAAAGCTCGCAGCAGCCGGAAGTCAGGCTCTGCCAGCTGAACTACGCGGACATAAAGGATTA
CTTGGCTCAGCCTGAGAAGATCACCTTGGTGTTCCTTGGAGTCGAGCTTGAATGAGAAAAGGCTCGCAC
GCCATGCCGGCGGAGTCCAGAGGGAGAGGAGGATGGTTTGGTGCATGGTTTGTCTCCGCATAGAAC
CTGGTGCAGCTGAGGAGTTAAGCAAAGACACGAAAAGTACTTTCTCCACCCGCAATGCCCGCTCT
TCTGCAGTTGAAAGAAAAGAGGCTGGGGTCTGGCTCAAGCAAGATCGGTGCTTGGCTGACAGTCTCGA
TACAAGTTCTGCCAACCTGTGGCAGCACAATAAAATTGAAGAAGGTGGCTACAAAAGAGTGTGTGTT
GAGAAAAGTGCCTAGTCTCCATGGCGTCCACAACACATCTTACCCGAGAGTTGATCCAGTCTGTAATTAT
GCAGGTCAATCATCCGGATGGAACCAAGTGTCTTTTAGGCCGGCAAAAAGAGATTTCCCCCGGCATGTT
ACTTGTCTCGTGGATTTATCGAGCCTGGGAGACAATAAGACGCTGTGCGGAGAGAGGTAGGAGAGG
AAAGTGGAGTCAAAGTTGGCCATGTTTCAAGTGTCTTGTGTCAGCCGTGGCAATGCCCTCTCTTAAT
GATTGGCTGTTTAGCTGTGGCAGTGTCTACAGAAATTAAGTTGACAAGAATGAAATAGAGGATGCCCGA
TGTTTCACTAGAGAACAGGTGGTGGATGTTCTTACCAAAGGGAAGCAGCAAGCATTTCTCGTCCACCAA
GCCGAGCTATTGCACATCAGTTAATCAAGCACTGGGTTGGCATGAACCCAGTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA




Protein Sequence: >RR201521 representing NM_001109010
 Red=Cloning site Green=Tags(s)

MSSVKRNPKREIISELHSSAAEGDVAKLGLSHSPSLLNETSENGWTALMYAARNGHPDAVQFLLEKGC
 DRSSVVKSRQTALDIAVFWGYKHIANLLANAKGGKPPWFLTSDVDGCENYFSRTLLDRRSDKRTNSDWLQ
 AKESHPTTVYILFSDLNPLVTLGGNKESQQPEVRLCQLNYADIKDYLAQPEKITLVFLGVELEMRKGS
 AHAGGVPEGEEDGLVAWFVLGIEPGAEEFKQRHENCYFLHPPMPALLQLKEKEAGVVAQARSVLAWHSR
 YKFCPTCGSTTKIEEGGYKRVCVRENCPSLHGVHNTSYPRVDPVVMQVIHPDGTCLLGRQKRFPPGMF
 TCLAGFIEPGETIEDAVRREVGEESGVKGVHVQYVSCQPWPMPSSLMIGCLAVAVSTEIKVDKNEIEDAR
 WFTREQVVDVLTGKGQQAFFVPPSRAIAHQLIKHWGMNPSL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001109010
ORF Size:	1386 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109010.1 , NP_001102480.1
RefSeq Size:	3880 bp
RefSeq ORF:	1389 bp
Locus ID:	367323
Cytogenetics:	9q36
MW:	51.4 kDa