

## Product datasheet for **RR201525L4V**

### **Nfil3 (NM\_053727) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Nfil3 (NM_053727) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nfil3
Synonyms:	E4BP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053727
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201525).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053727.2</a> , <a href="#">NP_446179.2</a>
RefSeq Size:	1990 bp
RefSeq ORF:	1389 bp
Locus ID:	114519
UniProt ID:	<a href="#">Q6IMZ0</a>
Cytogenetics:	17p14



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a transcriptional regulator that binds as a homodimer to activating transcription factor (ATF) sites in many cellular and viral promoters. The encoded protein represses Per1 and Per2 expression and therefore plays a role in the regulation of circadian rhythm. [provided by RefSeq, Feb 2014]