

Product datasheet for **RR201545**

Il17rb (NM_001107290) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Il17rb (NM_001107290) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Il17rb
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR201545 representing NM_001107290 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACGAGGACAGCCCTCTTTGTCTGTGAACCTCACCTCGCCAGGCTGCCTAAACCATGTAATGAAAT
ATAAGAAGCAGTGCATTGAAGCAGGAAGCCTGTGGGACCCAAACATCACAGCTTGAAGAAGAATGAGAA
GACCGTTGAAGTGAATTTACAACCAATCCCTTGGAAACAGATACATGGTTCTCATTTCGACGGGACACA
ATGTTGGGGTTTCTATAGTGTGGAGAACAACACTGACGAGGACATCTGTAGTCATTCCGGTCAATGATG
AGAGTGAAGGTGCCTTGTTGAGCTGACCCATATTTACACACCTGCGACAATGACTGCATTTCGACGCAA
AGGGACAGTGGTCTTTGCTCAGAGACAAGTGTCCCTCCCTCCAGATGACAACAGAAGCATGCTGAGA
GGCTGGCTGCCTCTCCTCCTGGTGTGCTGGTGGCTACATGGGTGCTGGCAGTTGGGATCTACCTAACTT
GGAGGCAAGGAAGGAGCACAAGACATCTTTCCATCACACAGTCTGTCGCTTACTGACTTTCTCAAACCTAC
TGCAGAAGTGAGTGCCTTGAATAATGCGAGAAAAGAAAATCGCCGAGATGGGGCCCGTGCAGTGGC
TTACCACTCAGAAGCAAGCAGCAGATAAAGTCGTCTTCTTCCCAGTGACGCCCTTCCCTTTGTGA
CAGTGCCTGTGGCCACAAGGAGGCAAGTCCACGAGAAATCTCAGGATCTGTTCCCTCTTGCCCTTAAC
CTCTTTGTAGCGATTTACGAGCCAGACTCATCTGCACAAATACCTGGTGGTCTATCTTGGGGGAGCAG
ACCTCAAAGGGGACTATAATGCCCTGAGAGTCTGCCCCAGTATCATCTCATGAAGGACGCCCCAGCTTT
CCACACAGAACTCTCAAGGCTACGCAAGCATGCCGCTGAAGAAACGGCCACAGGCCTGCCACGGCAGC
TGTTACCCCTTG

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR201545 representing NM_001107290
 Red=Cloning site Green=Tags(s)

MNEDSPSLSVNFTSPGCLNHVMKYKQKQIEAGSLWDPNITACKKNEKTVEVNFNTNSLGNRYMVLIRRD
 MLGVSIVLENKLRTRTSVIPVNDSEALVELTPYLHTCDNDCIRRKGTVVLCSETSPAPFPDDNRSMLR
 GWLPLLLVLLVATWVLAAGIYL TWRQGRSTKTSFPI TAMLLPLVKVLVVPSEICFHHTVCRFTDFLQNY
 CRSEVILEKWQKKKIAEMGPVQWLTTQKQAADKVVVFLPSDAPSLCDSACGHKEGSATENSQDLFPLAFN
 LFCSDFSSTHLHKYLVVYLGADLKGDYNALRVCPQYHLMKDAPAFHTELLKATQSMPLKKRPQACHGS
 CSPL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001107290

ORF Size:	1062 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107290.1 , NP_001100760.1
RefSeq Size:	2056 bp
RefSeq ORF:	1065 bp
Locus ID:	306247
Cytogenetics:	16p16
MW:	39.7 kDa