

Product datasheet for **RR201565L3V**

Adm (NM_012715) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adm (NM_012715) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Adm
Synonyms:	Ap; H39316; RATAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012715
ORF Size:	555 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201565).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012715.1 , NP_036847.1
RefSeq Size:	1395 bp
RefSeq ORF:	558 bp
Locus ID:	25026
UniProt ID:	P43145
Cytogenetics:	1q33
Gene Summary:	involved in the regulation of adrenal development; plays a role in inhibition of apoptosis [RGD, Feb 2006]



[View online »](#)