

Product datasheet for RR201580L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Bdkrb1 (NM_030851) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Bdkrb1 (NM_030851) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Bdkrb1

Synonyms: b1bkr; Bdkrb; BKR

Mammalian Cell

. . . .

Puromycin

Selection: Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_030851 **ORF Size:** 1011 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RR201580).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 030851.1, NP 110478.1

 RefSeq Size:
 1312 bp

 RefSeq ORF:
 1014 bp

 Locus ID:
 81509

 UniProt ID:
 P97583

Cytogenetics: 6q32

Gene Summary: may play a role in regulation of inflammatory response [RGD, Feb 2006]