

## Product datasheet for **RR201583L4V**

### Slc39a9 (NM\_001034929) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc39a9 (NM_001034929) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc39a9
Synonyms:	MGC125246
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001034929
ORF Size:	861 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201583).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001034929.1</a> , <a href="#">NP_001030101.1</a>
RefSeq Size:	1538 bp
RefSeq ORF:	864 bp
Locus ID:	314275
UniProt ID:	<a href="#">Q3KR82</a>
Cytogenetics:	6q24
Gene Summary:	May act as a zinc-influx transporter.[UniProtKB/Swiss-Prot Function]



[View online »](#)