

## Product datasheet for **RR201666**

### **Nt5c3a (NM\_001107862) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nt5c3a (NM\_001107862) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Nt5c3a  
**Synonyms:** Nt5c3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR201666 representing NM\_001107862  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGACTAATCAGGAGTCTGCTGTCGGTTTGAAATGATGCCAGAATTCAGAAAAGCTCAGTTCGTATCA  
AGAACCCAACGAGAGTAGAAGAAATTATCTGTGGTCTCATTAAAGGAGGCGCTGCCAACTTCAGATAAT  
AACCGACTTTGATATGACTCTAAGTAGGTTCTCCTACAATGGGAAAAGATGCCAACATGTCATAATATC  
ATTGACAACGTAAAGCTCGTTACAGATGAATGTCGAAGAAAGTTACTGCAGCTAAAGGAGCAGTACTACG  
CTATTGAGGTTGATCCTGTTCTCACCGTGAAGAGAAGTTCCTTACATGGTAGAATGGTATACCAAATC  
GCATGGTTTGCTTATTGAACAAGGCATACCAAAAGCTAAACTTAAAGAGATTGTGGCTGACTCTGATGTC  
ATGCTCAAGGAAGGATATGAGAAATTTCTTCGGTAACTTCAACAGCAGCGGTATCCCTGTGTTTATTTT  
CGGCTGGTATCGGTGACGACTTGAAGAAGTTATCCGGCAGGCTGGAGTTTACTCAACAGTCAAAGT  
AGTGTCTAACTTCATGGATTTTGTGAAAATGGAGTGTGAAAGGATTCAAAGGAGAGCTAATTCATGTA  
TTCAACAAACATGATGGTGCCCTGAAGAATACAGACTACTTCAGCCAACGAAAGACAACAGCAACATCA  
TTCTGCTGGGGGACTCCAAGGAGACTTAAGGATGGCCGATGGTGTAGCCAACGTCGAACACATTCTGAA  
AATTGGTTATCTGAATGACAGAGTAGACGAGCTTTTAGAAAAATACATGGACTCTTATGATATTGTTTTA  
GTAAAAGAAGAATCACTGGAGGTTGTCAACTCTATTTTACAGAAGACATTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR201666 representing NM\_001107862  
 Red=Cloning site Green=Tags(s)

MTNQESAVRLEMMPEFQKSSVRIKNPTRVEEIIICGLIKGGAAKLQIITDFDMTLSRFSYNGKRCPTCHNI  
 IDNCKLVTDECRRLQLKEQYYAIEVDPVLTVEEKFPYMVEWYTKSHGLLIEQGIPKAKLKEIVADSDV  
 MLKEGYENFFGKLQOHGIPVFI SAGIGDVL EEVIRQAGVYHSNVKVVSNFMDFDENGVLKGFKGELIHV  
 FNKHDGALKNTDYFSQLKDNSNIILLGDSQGLRMADGVANVEHILKIGYLNDRVDELLEKYMDSYDIVL  
 VKESLEVVNSILQKTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001107862

ORF Size: 891 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001107862.1</a> , <a href="#">NP_001101332.1</a>
<b>RefSeq Size:</b>	1617 bp
<b>RefSeq ORF:</b>	894 bp
<b>Locus ID:</b>	312373
<b>Cytogenetics:</b>	4q24
<b>MW:</b>	33.8 kDa