

Product datasheet for **RR201675L3V**

Gsn (NM_001004080) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gsn (NM_001004080) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gsn
Synonyms:	MGC95032
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001004080
ORF Size:	2340 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201675).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001004080.1 , NP_001004080.1
RefSeq Size:	2612 bp
RefSeq ORF:	2343 bp
Locus ID:	296654
UniProt ID:	Q68FP1
Cytogenetics:	3p11
Gene Summary:	plays a role in phosphoinositide mediated disassembly of actin filaments in Sertoli cell adhesion complexes; may regulate sperm release and turnover of the blood-testis barrier [RGD, Feb 2006]



[View online »](#)