

Product datasheet for **RR201698**

Setd7 (NM_001109558) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Setd7 (NM_001109558) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Setd7
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR201698 representing NM_001109558 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGATAGCGACGACGAGGTAGTGGAGGAAGCCGTGGAAGGGCACCTGGACGATGACGGATTACCACACG
GGTCTGCACTGTCACCTACTCCTCCACAGACAGATTTGAGGGCAACTTCGTGCATGGAGAAAAGAATGG
ACGTGGGAAATTCTTCTTTGACGGCAGCACCTGGAGGGATACTATGTGGATGATGCCCTGCAGGGC
CAGGGCGTTTACACCTATGAGGATGGAGGTGTTCTGCAAGGCACCTATGTGGATGGGGAGCTGAATGGGC
CAGCCCAGGAGTACGACACCGACGGGAGGCTGATCTTCAAGGGACAGTACAAAGACAACATTCGACACGG
AGTGTGTTGGATCCATTACCCAGATGGAGGCAGTCTCGTCGGGAAGTCAATGAAGATGGGGAGATGACA
GGGAAAAGATAGCCTACGTGTACCCTGATCAGAGGACTGCGCTCTATGGGAAGTTCAATGATGGAGAGA
TGCTGGAAGGGAAGCTGGCCACCCTCATGGCCACAGAAGAAGGAAGCCACACTTTGAAGTACATCGGG
AAGTCCGTGTTCCACTTTGATAAATCGACTTCCTCCTGCATCTCTAGCGATGCCCTCCTCCAGACCCT
TATGAGTCCGAGAGGGTTTATGTTGCCGACTCTCTCATCTCCAGTGCCGGAAGGACTGTTTTCAAGG
TAGCAGTGGGACCAATACTGTGATGTCATTTTATAATGGAGTCCGAATTACACCAAGAGGTTGACAG
CAGGGACTGGGCCCTAAATGGGAACACTCTCCTCCTGGATGAGGAAACAGTCATTGATGTGCTGAGCCC
TACAACCACGTGTCCAAGTACTGTGCCTCCTTGGGACACAAGGCGAATCACTCCTTTACTCCAACTGCG
TCTATGACATGTTTGTCCATCCCGTTTTGGGCCCATCAAGTGCATCCGCACCCTTCGAGCTGTGGAGGC
TGAGGAGGAGCTGACCGTTGCCTACGGCTATGACCACAGCCCCCGGGGAAGAGTGGGCCGAAGCCCC
GAGTGGTATCAGGTGGAGCTGAAGGCCTTCCAGGCCACCCAGCAAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR201698 representing NM_001109558
Red=Cloning site Green=Tags(s)


MDSDDVEVEEAEGHLLDDGLPHGFCTVYTSSTRDFEGNFVHGEKNGRGKFFFFDGGSTLEGYYVDDALQG
 QGVYTYEDGGVLQGTVDGELNGPAQEYDTDGRILFKGQYKDNIRHGVCWIHYPDGGSVLGEVNEDGEMT
 GEKIAYVYPDQRTALYKGFIDGEMLEGKLATLMATEEGKPHFELTSGSSVFHFDKSTSSCISDALLPDP
 YESERVYVADSLISSAGEGLFSKVAVGPNVTMSFYNGVRITHQEVDSDRWALNGNTLSLDEETVIDVPEP
 YNHVSKYCASLGHKANHSFTPNCVYDMFVHPRFGPIKCIRTLRAVEAEELTVAYGYDHSPPGKSGPEAP
 EYWOVELKAFQATOOK

TRTRPLEOKLISEEDLAANDILDYKDDDDKV


Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001109558

ORF Size:	1098 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001109558.1 , NP_001103028.1
RefSeq Size:	1977 bp
RefSeq ORF:	1101 bp
Locus ID:	689954
Cytogenetics:	2q26
MW:	40.5 kDa