

Product datasheet for **RR201722L3V**

Skap1 (NM_173311) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Skap1 (NM_173311) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Skap1
Synonyms:	Scap1; Skap55
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173311
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201722).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173311.2 , NP_775433.2
RefSeq Size:	1551 bp
RefSeq ORF:	1065 bp
Locus ID:	286975
UniProt ID:	Q4V7G1
Cytogenetics:	10q26-q31
Gene Summary:	cytoplasmic phosphoprotein; member of the src family of kinase enzymes; may associate with mast cell immunoreceptor signal transducer (MIST) [RGD, Feb 2006]


[View online »](#)