

Product datasheet for **RR201734L3V**

Fkbp1 (NM_001002818) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fkbp1 (NM_001002818) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fkbp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001002818
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201734).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001002818.2 , NP_001002818.1
RefSeq Size:	1300 bp
RefSeq ORF:	1044 bp
Locus ID:	406168
UniProt ID:	Q6MG81
Cytogenetics:	20p12
Gene Summary:	May be involved in response to X-ray. Regulates p21 protein stability by binding to Hsp90 and p21.[UniProtKB/Swiss-Prot Function]


[View online »](#)