

Product datasheet for **RR201740**

Acox2 (NM_145770) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Acox2 (NM_145770) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Acox2
Synonyms:	THCCox
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR201740 representing NM_145770
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGAGCCCAATGCACCGAGTGTCCCTGGGGACACCTGGAGCTGGCAAGTGCACCCGGACATAGACA
GTGAAAGGCACTCACCGTCCTTCAAGTGTGGAGCGACTCACCAACATCCTTGATGGAGGCCTCCCAAACAC
CGTGCTGCGAAGAAAAGTCGAAAGCATCATACAAAGTGACCCAGTGTTTAAATTTGAAGAAGCTTTACTTC
ATGACCCGAGAGGAGCTATATGAGGATGCGATTCAAAGAGATTCCATCTCGAGAAGCTAGCCTGGAGCC
TGGGCTGGTCAGAAGATGGTCTGAACGCATTTATGCTAACAGAGTCTTGATGGAAACGTCACCTTAAG
CTTACATGGTGTGCCATGAATGCTATCCGAAGCCTGGGCTCAGATGAACAGATTGCTAAATGGGGCCAA
CTCTGCAAAAACCTCAAATCATCACACATACGCCAGACAGAGCTGGGACACGGGACATACCTACAGG
GCCTGGAGACTGAAGCCACCTATGATGAAGCCAGGCAGGAGTTTGTGATACACAGCCCTACGATGACTTC
CACCAAGTGGTGGCCTGGGGACTTGGGATGGTCCGTCACCCATGCTGTGGTCTAGCCAGTTGACCTGC
TTAGGAGTCCGGCACGGCATGCACGCCTTCATTGTGCCATTTCGAGCCTAGAGGATCACACCCCACTGC
CAGGAATCACAGTTGGGGACATAGGCCCAAGATGGGTTTGAACACATAGACAATGGCTTCTGCAACT
GAACCAGTGGGGTTCCAGAGAAAACATGCTCAGTCGCTTTCAGAGGTCTTGCCAGATGGTACCTAC
CAGAGGCTTGGGACGCCAGAGCAATTATCTTGGCATGTTGGTGACCCGGTGCAGCTGCTGTGTAAG
GAATCCTACCCTCCCTCCAGAAGCTTGTCATATTGCCACGCGCTACTCAGTAATCCGCCATCAGTCTCG
ACTTCGGCCAGTGACCCAGAGGCAAAAATCCTGGAATACCAGACGCAGCAGCAGAAAACCTTCTCCTCAG
CTTGCTGTGAGCTATGCCTTCCACTTCAGGCCACCAGCCTCTCAGAATCTTCCACAGCTCCTACAGTG
CTATTCTGAAGAGAGACTTCAGCCTCTGCCTGAGCTCCATGCATTGAGCACTGGTATGAAGGCCACGTT
TGCAGACTTCTGTGCCAGGGGGCCGAGATCTGTGCGAGAGCTTGTGGGGCCATGGCTACTCAAAGCTG
AGCGCCCTGCCGACACTGGTTGCTCGAGCAACAGCCTCTTGCACATATGAGGGTGAGAATACGGTGTCT
ACCTGCAAGTGGCCAGGTTTCTGATGAAGAGCTATCTGCAGGCTCAAGCGTCCCCAGGCGCCACACCACA
GAAGCCTCTCCCTCAGTCCGTATGTATATTGCCACACAAAGGCCAGCCAGGTGCTCAGCCAGACTGCA
GCTGACTTCCGCTGCCAGATGTCTATACCACAGCCTGGGCATATGTGTCTACCAGGCTCATAAGAGATG
CAGCACACCGTACACAGACCCTCATGAAGTCCGGGTTGACCAGCATGATGCCTGGAATCAGACTACTGT
CATCCACCTTCAGGCTGCTAAGGCTTACTGCTACTTCATCACTGTGAAGAATTTCAAGGAAGCTGTGGAG
AAACTAGACAAGGAACCAGAGATTCAGCGTGTGCTCCAACGCCTCTGTGACCTCTATGCCTTACACGGTG
TTCTGACTAACTCAGGGACTTTCTGCATGATGGCTTCTGTCTGGGGCCAGGTGGACATGGCCAGAGA
AGCCTTCTAGACCTGCTTCCCTTGATCCGGAAGGATGCCATCTTGTTAACCGATGCTTTTGACTTCTCG
GACCATTGTTTAAACTCGGCACTTGGCTGTTATGATGGACACGTCTACGAACGCCTGTTTGAGTGGGCTC
AGAAGTCCCAGCCAATACTCAGGAGAACCCTGCCTATAAGAAGTATATCCGACCACTGATGCTCGGCTG
GAGACACAAGATG

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RR201740 representing NM_145770
 Red=Cloning site Green=Tags(s)


MGSPMHRVSLGDTWSQVHPDIDSERHSPSF SVERL TNILDGGLPNTVLRKVESIIQSDPVFNLKKLYF
 MTR EEL YEDAIQKRFHLEKLAWSL G WSE D G P E R I Y A N R V L D G N V N L S L H G V A M N A I R S L G S D E Q I A K W G Q
 L C K N F Q I I T T Y A Q T E L G H G T Y L Q G L E T A Y D E A R Q E F V I H S P T M T S T K W W P G D L G W S V T H A V V L A Q L T C
 L G V R H G M H A F I V P I R S L E D H T P L P G I T V G D I G P K M G L E H I D N G F L Q L N H V R V P R E N M L S R F A E V L P D G T Y
 Q R L G T P Q S N Y L G M L V T R V Q L L C K G I L P S L Q K A C I I A T R Y S V I R H Q S R L R P S D P E A K I L E Y Q T Q Q K L L P Q
 L A V S Y A F H F T A T S L S E F F H S S Y S A I L K R D F S L L P E L H A L S T G M K A T F A D F C A Q G A E I C R R A C G G H G Y S K L
 S G L P T L V A R A T A S C T Y E G E N T V L Y L Q V A R F L M K S Y L Q A Q A S P G A T P Q K P L P Q S V M Y I A T Q R P A R C S A Q T A
 A D F R C P D V Y T T A W A Y V S T R L I R D A A H R T Q T L M K S G V D Q H D A W N Q T T V I H L Q A A K A Y C Y F I T V K N F K E A V E
 K L D K E P E I Q R V L Q R L C D L Y A L H G V L T N S G D F L H D G F L S G A Q V D M A R E A F L D L L P L I R K D A I L L T D A F D F S
 D H C L N S A L G C Y D G H V Y E R L F E W A Q K S P A N T Q E N P A Y K K Y I R P L M L G W R H K M

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Cloning Scheme:



Plasmid Map:


ACCN: NM_145770

ORF Size: 2043 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145770.2](#), [NP_665713.2](#)

RefSeq Size: 2504 bp

RefSeq ORF: 2046 bp
Locus ID: 252898
Cytogenetics: 15p14
MW: 76.7 kDa
Gene Summary: plays a role in fatty acid metabolism [RGD, Feb 2006]