

Product datasheet for **RR201743**

Fshb (NM_001007597) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fshb (NM_001007597) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fshb
Synonyms: FSH-B; FSH-beta
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201743 representing NM_001007597
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGATGAAGTCGATCCAGCTTTGCATCCTACTCTGGTGCTTGAGAGCAGTCTGCTGCCATAGCTGTGAAC
TGACCAACATCACCATCTCAGTAGAGAAGGAAGAGTGCCGTTTCTGCATAAGCATCAATACCACTTGGTG
TGAGGGCTACTGCTACACCAGGGATCTGGTGTATAAGGACCCAGCTAGACCAAACACCCAGAAAGTATGC
ACTTCAAGGAGCTGGTGTACGAGACCATAAGATTGCCTGGCTGTGCCCGCCACTCAGACTCCCTCTACA
CATATCCAGTAGCCACTGAATGCCACTGTGGCAAGTGTGATAGTGACAGCACCGACTGCACTGTAAGAGG
CCTGGGACCCAGCTACTGCTCCTTCGGTAAATGAAAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR201743 representing NM_001007597
Red=Cloning site Green=Tags(s)

MMKSIQLCILLWCLRAVCCHSCELTNITISVEKEECRFCISINTTWCEGYCYTRDLVYKDPARPNTQKVC
TFKELVYETIRLPGCARHSDSLYTPVATECHCGKCDSDSTDCTVRGLGPSYCSFGEMKE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001007597

ORF Size: 390 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001007597.2 , NP_001007598.1
RefSeq Size:	393 bp
RefSeq ORF:	393 bp
Locus ID:	25447
UniProt ID:	P18427
Cytogenetics:	3q33
MW:	14.8 kDa
Gene Summary:	stimulates development of follicle and spermatogenesis in the reproductive organs [RGD, Feb 2006]