

Product datasheet for **RR201753L4V**

Zmat3 (NM_022548) Rat Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Zmat3 (NM_022548) Rat Tagged ORF Clone Lentiviral Particle |
| Symbol: | Zmat3 |
| Synonyms: | PAG608; Wig1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_022548 |
| ORF Size: | 867 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RR201753). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_022548.1 , NP_071993.1 |
| RefSeq Size: | 7513 bp |
| RefSeq ORF: | 870 bp |
| Locus ID: | 64394 |
| UniProt ID: | O08781 |
| Cytogenetics: | 2q24 |
| Gene Summary: | nuclear zinc finger protein; may play a role in facilitating the apoptotic activities of p53 [RGD, Feb 2006] |



[View online »](#)