

Product datasheet for **RR201757**

Rab26 (NM_133580) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rab26 (NM_133580) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rab26
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR201757 representing NM_133580
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGCTGGTGGGGATTCCGGTGTGGGAAGACGTGCCTGCTCGTGCCTTCAAGGATGGGGCCTTCCTGG
 CTGGTACCTTCATCTCCACTGTGGGCATCGACTTCCGGAATAAAGTTCTGGATGTGGATGGCATGAAGGT
 GAAGCTGCAGATCTGGGACACGGCTGGTCAGGAGCGGTTCCGAAGTGTCACCCATGCCTACTACCGGGAT
 GCTCACGCACTGTTGCTGCTCTACGACATCACCAACAAAGATTCCCTCGACAACATCCAGGCCTGGTTGA
 CAGAAATCCAGGAATATGCCAGCAGGACGTGGTCTCATGCTGTTGGGGAACAAGTTGACTCTACTCA
 AGAACGAGTGGTAAAGAGGGAAGATGGGGAGAAATTAGCCAAGGAGTATGGGCTGCCATTCATGGAGACC
 AGCGCCAAGAGCGGCCTCAATGTGGACTTGGCTTTTACAGCCATAGCAAAGGAGCTGAAACAAAGATCCA
 CCAAGGCTCCCAGTGAGCCGCGCTTCAGGCTGCATGACTATGTGAAGAGGGAGGGCCGAGGGTCTCCTG
 CTGTCGACTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR201757 representing NM_133580
Red=Cloning site Green=Tags(s)

MLVGDGSGVGTCLLVRFKDGAFLAGTFISTVGIIDFRNKVLDVDGMKVKLQIWDTAGQERFRSVTHAYYRD
 AHALLLLYDITNKDSFDNIQAWL TEIQEY AQQDVVLM L L GNKVDSTQERVVKREDGEKLAKEYGLPFMET
 SAKSGLNVDLAFATAIAKELKQRSTKAPSEPRFRLHDYVKREGRGVSCRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_133580

ORF Size: 570 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133580.1</u> , <u>NP_598264.1</u>
RefSeq Size:	1098 bp
RefSeq ORF:	573 bp
Locus ID:	171111
UniProt ID:	<u>P51156</u>
Cytogenetics:	10q12
MW:	21.5 kDa
Gene Summary:	may play a role in vesicular transport [RGD, Feb 2006]