

## Product datasheet for **RR201852**

### Cyb5b (NM\_030586) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cyb5b (NM\_030586) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cyb5b  
**Synonyms:** Cyb5m; omb5  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR201852 representing NM\_030586  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCGACCCAGAGGCCAGCGGCAGCGGAGGAATGGCCAGGGCTCCGACCCCGCCGTCACTACTACC  
GGCTGGAGGAGGTGGCGAAGCGAAATACCGCGGAGGAGACGTGGATGGTATCCACGGGCGAGTCTACGA  
TATCACCCGCTTCTCAGCGAGCATCCTGGTGGAGAAGAGTTCTACTGGAACAAGCTGGTGCCGATGCA  
ACTGAAAGCTTTGAAGATGTCGGCCACTCCCCTGATGCCAGGGAGATGCTAAAGCAGTACTACATTTGGT  
ATGTCCATCCGAATGACCTTAAACCTAAAGATGGTGACAAGGACCTTCAAAAAACAATTCATGCCAAAG  
TTCTGGGCATATTGGATTGTCCTCATCGTGGGTGCTATTCTTATAGTTTCTGTATCGTCACTTCTGG  
GCTGACAGCAAATCCTCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR201852 representing NM\_030586  
**Red=Cloning site Green=Tags(s)**

MATPEASGSRNGQSDPAVTYYRLEEVAKRNTAEETWMVIHGRVYDITRFLSEHPGGEEVLLAQAGADA  
TESFEDVGHSPDAREMLKQYYIGDVHPNDLKPDKGDKDPSKNNSCQSSWAYWIVPVGAILIGFLYRHFV  
ADSKSS

**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_030586

ORF Size: 438 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_030586.2</a></u> , <u><a href="#">NP_085075.1</a></u>
<b>RefSeq Size:</b>	2745 bp
<b>RefSeq ORF:</b>	441 bp
<b>Locus ID:</b>	80773
<b>UniProt ID:</b>	<u><a href="#">P04166</a></u>
<b>Cytogenetics:</b>	19q12
<b>MW:</b>	16.3 kDa
<b>Gene Summary:</b>	a cytochrome with isoforms located in endoplasmic reticulum and mitochondrial outer membranes [RGD, Feb 2006]