

Product datasheet for **RR201867**

Ifit2 (NM_001024753) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ifit2 (NM_001024753) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ifit2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR201867 representing NM_001024753
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGTACAGCCAGTAAGGAATCATTGGAGAGCAGTCTACAGCAGCTCAAATGCCACTTCACCTGGAACC
 TGAGAGCAGAAGACGAGTCCCTTGGATGAGTTTGAGGACAGGGTGTAAACAAAGATGAGTTTCAGAACAG
 TGAGTTTAAAGCCACCATGTGCAACCTACTGGCCTATGTTAAGCATTGCAGAGGTCTAAATGAGGCAGCA
 CTGAAATGCTTAGGGGAAGCCGAAGACTTCATCCGGCAACAACATCCTGATCAAATAGAAATCAAAAGTC
 TGGTCACCTGGGAAACTATGCTTGGGTTTACTATCACATGGGCCAGCTGTCCAAAGCTCAGGAGTATCT
 GGACAAGGTGAACAGGTCTGTAAGAAGTTCTCCAGTCCCTACAGGATTGAGAATCCTGTGCTTGACTGT
 GAGGAGGGGTGGCCCGACTAAGTGTACCAAGAACCAAAATGAGAGAATGAAGGTGTGCTTCGAGAAGG
 CTCTGGAAGGACCCAAAGAACCAGAATCCACTTCTGGATGGCCATTGCAAACCTACCGTCTGGATGA
 CTGGCCAGCATCCAATGACTATATTGACTCTCTGGAGCAAGCCATTAGCCTGTCTCTGACAACACGTAT
 GTTAAAGTCTCTTGGCAATGAAACTTGAGGAGGTGCATGAAAACCGAGCAAAGGAGCTAGTTGAAGAAG
 CCTTGAAGAAGGACCCAAGTGAATAGATAACAATGCTTGGTGCAGCCAAGTTTTACGTCAAGGTACATGA
 CACAGACAGAGCAATACAGTTGCTTAAAAAGGCTTTGGAATCCATGCCCAACAATGCTTATGTGCATTAC
 TACCTTGGGTGCTGTTACAGGTCAAAAGTCTTACATACTGAACACCGAAGAAACCACATCTAATGGAA
 ATAGGGAGAAGTTAGAGGCAGTACATAACCAGCACTGAATCACTTACGGAAAGCAGAGGATATCAAGGA
 GATGATCGAAAATTCCTGTAGTCACCTCGCCGGCCTTTATGTACAACAAAGCAGTATGAAGAAGCCGAT
 TATTACTCCAGAAAGAATTCAACAAGGACCTCCGTCCTTAAACAGTTACTCCATCTACGGTATG
 GCAACTTTCAGTTCTACGAAAGGAAGTGAAGACAGAGCCATCTACCACTATATGGAGGGTGTGAAAAT
 AAATCGGGAATCAAAGCCTAAAGGAAAGATGAGAGACAAACTTCGAAAACCTTGCTCTAAGCGCGCTTCT
 AGAGATGAATCTGATCCTGAGGCCTTGCATCTTGGCATTCTTAGGGAGCATAGAGAAGGTGAGGGAG
 CTGATAAAGCCTCGGAGGGAGAAGAGGACCCTGGAAACCGGGTCCCTTACGATCTTTGAACTGATAAA
 GCCTCGGAGGGAGAAGAGGACCCTGAAACCGGGTCCCTTACGATCTTTGGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR201867 representing NM_001024753
 Red=Cloning site Green=Tags(s)


MSTASKESLESSLQQLKCHFTWNLRAEDESLEDEFEDRVFNKDEFQNSEFKATMCNLLAYVKHCRGLNEAA
 LKCLGEAEDFIRQQHPDQIEIKSLVTWGNyawvyyHMQLSKAQEYLDKVKQVCKKFSSPYRIENPVLDC
 EEGWARLKCTKNQNERMKVCFEKALEKDPKNPESTSGWAIANYRLDDWPASNDYIDSLEQAISLSPDNTY
 VKVLLAMKLEEVHENRAKELVEEALKKDPSAIDTMLGAAKFYVKVHDTDRAIQLLKKALESMPNNAYVHY
 YLGCCYRSKVLHILNTEETTSGNREKLEALIQPALNHLRKAEDIKEMIENSCSHLAGLYVTTKQYEEAD
 YYFQKEFNKDLRPLPKQLLHLRYGNFQFYERKCEDRAIYHYMEGVKINRESKPKGMRDKLRKLALRRLS
 RDESDPEALRILAFLEHREGQGADKASEGEEDPGRVPSASLELIKPREKRTLETGSLQHLW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

Plasmid Map:


ACCN: NM_001024753

ORF Size: 1392 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq Size: 3433 bp

RefSeq ORF: 1395 bp

Locus ID: 294091

Cytogenetics: 1q53

MW: 53.9 kDa