

Product datasheet for **RR201870**

Mup5 (NM_203325) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mup5 (NM_203325) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mup5
Synonyms:	PGCL6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR201870 representing NM_203325 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAAGCTGTTGCTGCTGCTGTGTCTGGCCTGACCCTGGTCTGTGGCCATGAAGAAGAAGCTAGTT
 TCGAGAGAGGGAACCTCGATGTGGACAAGCTCAATGGGATTGGTTTCTATTGTCATGGCCTCTGACAA
 AAGAGAAAAGATAGAAGAGAATGGCAGCATGAGAGTTTTGTGCAGCACATCGATGCTTGGAGAATTCC
 TTAGGCTTCAAGTTTCATATTAAGGTAATGGAAAGTGCAGGGAACCTCTATTTGGTTGCCTACAAAACGC
 CAAAGGATGGCGAATATTTTGTGAATATGACGGAGGGAATACATTTAATACTTAAGACAGACTATGA
 CAGATATGTCATGTTTCATCTCGTTAATTTCAAGGACGGGAAACCTTCCAGCTGATGAATCTCTTAGGC
 AGAACAAAGGATCTGAGTTCAGATATCAAGGAAAAGTTTGCAAACTATGTGTGGCGCATGGAATCACTA
 GGGAGAATATCATTGACGTAACCAAGACTGATCGCTGTCTCCAGGCCGAGGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RR201870 representing NM_203325 Red=Cloning site Green=Tags(s)
-------------------	--

MKLLLLLLCLGLTLVCGHEEEASFERNLDVDKLNQDWFISIVMASDKREKIEENGSMRVFVQHIDVLENS
 LGFKFHIKVNKCRELYLVAYKTPKDGEYFVEYDGGNTFNILKTDYDRYVMFHLVNFKDGETFQLMNLG
 RTKDLSSDIKEKFAKLCVAHGITRENIIDVTKDRCLQARG

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:	SgfI-MluI
--------------------	-----------


[View online »](#)

Cloning Scheme:



ACCN: NM_203325

ORF Size: 543 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_203325.1](#), [NP_976070.1](#)

RefSeq Size: 789 bp

RefSeq ORF: 546 bp


Locus ID: 298107

Cytogenetics: 5q24

MW: 20.8 kDa

Gene Summary: member of the alpha(2u)-globulin multigene family; may be involved in pheromone response [RGD, Feb 2006]

Product images:



Circular map for RR201870