

Product datasheet for **RR201879L4V**

Haus8 (NM_001024971) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Haus8 (NM_001024971) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Haus8
Synonyms:	Ny-sar-48; RGD1311565
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024971
ORF Size:	1155 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201879).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024971.2 , NP_001020142.2
RefSeq Size:	1464 bp
RefSeq ORF:	1158 bp
Locus ID:	290626
UniProt ID:	Q5BK57
Cytogenetics:	16p14
Gene Summary:	Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.[UniProtKB/Swiss-Prot Function]


[View online »](#)