

## Product datasheet for **RR201885L3V**

### Necap1 (NM\_001029919) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Necap1 (NM_001029919) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Necap1
Synonyms:	RGD1306053
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001029919
ORF Size:	831 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201885).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001029919.2</a> , <a href="#">NP_001025090.1</a>
RefSeq Size:	2459 bp
RefSeq ORF:	834 bp
Locus ID:	312694
UniProt ID:	<a href="#">P69682</a>
Cytogenetics:	4q42
Gene Summary:	Involved in endocytosis.[UniProtKB/Swiss-Prot Function]


[View online »](#)