

Product datasheet for **RR201889L4V**

Strn3 (NM_001029897) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Strn3 (NM_001029897) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Strn3
Synonyms:	Gs2na; Sg2na
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001029897
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201889).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001029897.2 , NP_001025068.1
RefSeq Size:	4152 bp
RefSeq ORF:	2385 bp
Locus ID:	114520
UniProt ID:	P58405
Cytogenetics:	6q23
Gene Summary:	human homolog is an S, G2 phase nuclear autoantigen detected in a patient with lung and bladder cancer [RGD, Feb 2006]


[View online »](#)