

Product datasheet for **RR201890**

Olr251 (NM_001000218) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr251 (NM_001000218) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olr251
Synonyms:	ratchr1-165674861-165673917_ORF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR201890 representing NM_001000218 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCTTTGGCTAGAGAATGGGAACCACTACAGTGGCAGAGTTCCTTTATTGGGATTAACAGATGACC
 CAGTTCCTTAGAGTCATTCTCTTACCATCATTTTGTGCATCTACCTGGTACTGTGTCTGGGAACCTCAG
 CACCATCCTTCTCATCAGAGTCTCCTCCAGCTCCATCACCCATGTACTTTTTCTCAGCCACTTGGCT
 TCTGTTGACATAGGCTACTCATCTTCTGTACACCCAATATGCTTGTTAATTTCTGGTAGAGAGAAATA
 TCATCTCTACCTGGGATGTGGCATAACAGCTCAGCCTAGTTGCTTTTTTGGTGCAGCTGAATGCTTCCT
 TCTGGCTGCCATGGCTTATGATCGCTTCATGGCAATCTGTAACCCACTGCTTTATTCTACCAAAATGTCC
 ACACAAGTCTGTATCCTGTTGGTTGTGGCTTCTTATATAGGGGGTTTCCTTAACGCTTCCTTTGCCATGG
 TTTCTTTTGTCTTTTCTTCTGTGGACCAATATAATCAATCATTTTTCTGTGATTTTGCTCCTTT
 GATTGAGCTCTCCTGTTCTGATGTCAAGTGTATGTTATTGCTACCTCATTTTCTGCTGGCTCAGTACT
 ATGATCACAGCACTTGTATAGCCATCTCCTATTCCTATATTCTCATTACCATCTGAAGATGCACTCCA
 CTGAGGGTCGACAGAAGGCCTTCTCCACCTGCACCTCCCATCTCACTGCAGTCACTTTATTCTATGGAAC
 TGTTACATTCATTTATGTGATGCCCAAGTCCAGTTTCTCCATTGACCAGAACAGATAGTGTCTGTCTTC
 TACTTGTGGTGATCCCATGTTGAACCCCTCATTTACAGTCTCAGAAATAATGAAATTAAGTGCTC
 TGAAGAAACAGATTTGGAGGAAACAGTTTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RR201890 representing NM_001000218
 Red=Cloning site Green=Tags(s)


MSWLENGNHTTVAEFLLLGLTDDPVLRLVILFTIILCIYLVTVSGNLSTILLIRVSSQLHHPMYFFLSHLA
 SVDIGYSSSVTPNMLVNFLVERNIISYLGCGIQLSLVAFFGAAECFLAAMAYDRFMAICNPLLSTKMMS
 TQVCILLVVASYIGGFLNASFAMVSFCSFLFCGPNIINHFFCDFAPLIELSCSDVSYYVIATSFSSAGSVT
 MITALVIAISYSYILITILKMHSTEGRQKAFSTCTSHLTAVTLFYGTVTFIYVMPKSSFSIDQNKIVSVF
 YLLVIPMLNPLIYSLRNNEIKSALKKQIWRKTVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001000218

ORF Size: 942 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min.
	2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
	3. Close the tube and incubate for 10 minutes at room temperature.
	4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
	5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000218.1, NP_001000218.1</u>
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	293394
Cytogenetics:	1q33
MW:	34.8 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RR201890