

Product datasheet for **RR201949**

Hfe (NM_053301) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hfe (NM_053301) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hfe
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR201949 representing NM_053301 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACCGATCAGCTGGGCTCCCTGTGCGGCTGCTATTGCTGCTGCTGTTGTTGCTGCTGTGGTCCGTTG
CCCCGACGGCGCTGCGGCCGGTTCACATTCTCTACGATATCTTTCATGGGTGCCTCAAAGCCAGACCT
CGGGCTGCCTTTCTTTGAGGCTCTGGGTTATGTGGATGACCAGCTCTTTGTATCCTACAATCAGGAGAT
CGCCGTGCTGAGCCAGGGCCCATGGATCTTGGGGCAGACATCAAGCCAGCTGTGGCTGCAGCTGAGTC
AGAGCCTGAAAGGGTGGGATTACATGTTTCATAGTGGACTTCTGGACCATCATGGGCAACTATAACCACAG
TAAGGTCACGAAGTTGAGAGTGGTGCCTGAGTCTCACATCCTGCAGGTGATCCTAGGATGTGAGGTGCAT
GAAGACAACAGTACCAGTGGCTTCTGGAAATATGGCTACGATGGGCAAGATCACCTTGAATTCTGCCCA
AGACACTGAACTGGAGTGCAGCCGAGCCAAGGGCCTGGGCCACCAAGATGGAGTGGGAAGAGCACAGGAT
CCGTGCCAGACAGAGCAGGGACTACCTGCAGAGGGACTGCCCCAGCAGCTGAAGCAGGTCTGGAGCTC
CAGAGAGGGTCTGGGACAGCAAGTGCCTACTTTGGTAAAAGTACTCGCCACTGGGCCTCTACAGGGA
CCTCTCTAAGGTGTCAGGCTCTGAATTTCTCCCCAGAACATCACTATGAGGTGTTGAAGGACAGCCA
GCCCTAGATGCAAGGATGTCAACCTGAGAAGTGTGCCAAATGGGGATGGGACCTATCAGGGCTGG
CTGACCTTGCTGTGGCCCTGGAGAAGAGACAAGTTTCAGCTGTCAAGTGGAGCACCCAGGCCTGGATC
AGCCTCTCACTGCCACTTGGGAGCCCTCAGGTCTCAGGACATGATTATTGGAATCATAAGTGGGATCAC
CATTTGTGCCATCTTCTTTGTTGGAATTCTGATCCTAGTCTTAAGGAAAAGGAAGTTTCAGGAGGAACC
ATGGGTGACTATGTCTTAACAGAGTGTGAGATGGGTGACTATGTCTTAACAGAGTGTGAG

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR201949 representing NM_053301
 Red=Cloning site Green=Tags(s)

MDRSAGLPVRLLLLLLLLLLWSVAPQALRPGSHSLRYLFMGASKPDLGLPFFFEALGYVDDQLFVSYNHES
 RRAEPRAPWILGQTSSQLWLQLSLSLKGWDYMFIVDFWTIMGNYNHKSVTKLRVVPESHILQVILGCEVH
 EDNSTSGFWKYGYDGDHLEFCPKTLNWSAAEPRAWATKMEWEHRIRARQSRDYLQRDCPQQLKQVLEL
 QRGVLGQQVPTLVKVTRHWASTGTSLRCQALNFFPQNI TMRWLKDSQPLDAKDVPENVL PNGDGT YQGW
 LTLAVAPGEETRFSCQVEHPGLDQPLTATWEPSRSQDMIIGIISGITICAIFFVIGILILVLRKRKVS
 GGTMGDYVLTCEMGDYVLTCE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053301

ORF Size:	1080 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	1670 bp
RefSeq ORF:	1083 bp
Locus ID:	29199
UniProt ID:	O35799
Cytogenetics:	17p11
MW:	41 kDa
Gene Summary:	membrane protein; may be involved with iron metabolism; defects in human homolog lead to hereditary hemochromatosis, an iron storage disorder [RGD, Feb 2006]