

Product datasheet for **RR201996L4V**

Serping1 (NM_199093) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serping1 (NM_199093) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Serping1
Synonyms:	MGC72585
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199093
ORF Size:	1512 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR201996).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199093.1 , NP_954524.1
RefSeq Size:	1850 bp
RefSeq ORF:	1515 bp
Locus ID:	295703
UniProt ID:	Q6P734
Cytogenetics:	3q24
Gene Summary:	human homolog acts as a serine protease inhibitor; involved in regulation of the classical complement pathway [RGD, Feb 2006]



[View online »](#)