

Product datasheet for **RR202020L4V**

Lifr (NM_031048) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lifr (NM_031048) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lifr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031048
ORF Size:	3279 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202020).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031048.1 , NP_112310.1
RefSeq Size:	3407 bp
RefSeq ORF:	3282 bp
Locus ID:	81680
UniProt ID:	O70535
Cytogenetics:	2q16
Gene Summary:	alpha chain of leukemia inhibitory factor glycoprotein and growth factor [RGD, Feb 2006]



[View online »](#)