

Product datasheet for **RR202023L4V**

St8sia2 (NM_057156) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St8sia2 (NM_057156) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	St8sia2
Synonyms:	Siat8b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_057156
ORF Size:	1125 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202023).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_057156.1 , NP_476497.1
RefSeq Size:	1128 bp
RefSeq ORF:	1128 bp
Locus ID:	117523
UniProt ID:	Q07977
Cytogenetics:	1q31
Gene Summary:	human homolog catalyzes the polysialylation of the neural cell adhesion molecule NCAM and may play a role in neurogenesis [RGD, Feb 2006]



[View online »](#)