

Product datasheet for **RR202041L3V**

Calhm2 (NM_001008306) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Calhm2 (NM_001008306) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Calhm2
Synonyms:	Fam26b; RGD1308276
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008306
ORF Size:	969 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202041).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008306.1 , NP_001008307.1
RefSeq Size:	1855 bp
RefSeq ORF:	972 bp
Locus ID:	294019
UniProt ID:	Q5R1Q8
Cytogenetics:	1q54
Gene Summary:	Pore-forming subunit of a voltage-gated ion channel.[UniProtKB/Swiss-Prot Function]



[View online »](#)