

## Product datasheet for **RR202046**

### Tspan13 (NM\_001013244) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tspan13 (NM\_001013244) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tspan13  
**Synonyms:** Tm4sf13  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR202046 representing NM\_001013244  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGTGTGCGGGGCTTCTCGTGCTCCAAGAACTGCCTGTGCGCGCTCAACCTGCTATACACCCTGGTCA  
 GTCTGCTGCTCATTGGAATTGCTGCTTGGGGCATCGGCTTTGGGCTGATCTCCAGTCTCCGAGTGGTGGG  
 CGTGGTCAATCGCTGTGGGCATCTCCTGTTCTGATTGCCTTGGTGGGGCTGATTGGAGCTGAAAAAC  
 CATCAGGTGTTGCTATTTTTATATGATCATTCTCTACTCGATTTATTGTCCAGTTTTCGGTATCTT  
 GTGCTTGTGTTAGCTCTGAATCGGGAGCAACAGGGTCAGCTCCTGGAGGTTGGATGGAACAATACAGCAAG  
 TGCTAGAAATGACATCCAGAGGAATTTGAACTGCTGTGGCTTCCGAAATTATAACCCAAACGACACCTGC  
 CCGGCTAGCTGTGCTAAAAGTAACCAGAAGTGCCTCTGTTGCCCCATAATCGGAGAGTATGCTGGAG  
 AGGTTTTGAGATTTGTTGGTGGCATTGGCCTCTTCTTCAAGTTTACAGAGATCCTGGGTGTTGGCTGAC  
 ATACAGATACCGGAACCAGAAAGACCCTCGGCCAACCCGAGCGCTTTCCTT

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR202046 representing NM\_001013244  
 Red=Cloning site Green=Tags(s)

MVCGGFSCSKNCLCALNLLYTLVSLLLIGIAAWGIGFGLISSLRVVGVIIVGIFLFLIALVGLIGAVKH  
 HQVLLFFYMIILLVFIQFVSVSACLALNREQQQLLEVGNNTASARNDIQRNLNCCGFRNYPNDTC  
 PASCASNQKCSSCAPIIGEYAGEVLRVGGIGLFFSFTEILGVWLTYYRNRQKDPANPSAFL

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:** SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM\_001013244

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001013244.1](#), [NP\\_001013262.1](#)

**RefSeq Size:** 1943 bp

**RefSeq ORF:** 615 bp

**Locus ID:** 366602

**UniProt ID:** [Q5FVL6](#)

**Cytogenetics:** 6q16

**MW:** 22.3 kDa