

Product datasheet for **RR202060L4V**

Trappc4 (NM_001003708) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trappc4 (NM_001003708) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trappc4
Synonyms:	Sbdn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001003708
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202060).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001003708.1
RefSeq Size:	863 bp
RefSeq ORF:	660 bp
Locus ID:	367073
UniProt ID:	Q69BT7
Cytogenetics:	8q22
Gene Summary:	mouse homolog binds syndecan-2 and may play a role in synaptic vesicle transport and dendritic spine formation [RGD, Feb 2006]



[View online »](#)