

Product datasheet for **RR202119**

Mvk (NM_031063) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Mvk (NM_031063) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Mvk
 Synonyms: Lrbp
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR202119 representing NM_031063
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTGTCAGAAGTCTGCTGGTGTCTGCTCCAGGAAAGTCATTCTCCATGGAGAACATGCTGTGGTCC
 ATGGCAAGGTAGCTCTGGCGGTGGCCTTGAAGTGAAGAACATTTCTCGTCTGCGACCGCAGAGCAATGG
 GAAAGTGAGCCTCAATTTACCAAACGTCGGTATTAAGCAGGTCTGGGATGTGGCCACACTTCAGCTGCTG
 GACACAGGCTTTCTTGTAGCAAGGCGATGTCCCGGCACCCACCTTGAGCAACTGGAGAAGCTGAAGAAGG
 TGGCGGGCCTCCCCGAGACTGTGTAGGCAACGAAGGCCTGTCTCTGCTTGCCTTCTGTACCTGTACCT
 GGCTATCTGCCGAAACAGAGGACACTCCCAAGCCTGGACATCATGGTGTGGTGGAACTGCCCCCTGGG
 GCGGGCTTGGGCTCCAGTGCAGCCTACTCGGTGTGTGTGGCAGCCGCCCTCCTGACTGCCTGTGAGGAGG
 TCACCAACCCGCTCAAGGACAGGGGCTCCATTGGCAGTTGGCCGAGGAGGACCTGAAGTCAATTAACA
 GTGGGCTACGAGGGGGAGAGAGTATCCATGGGAACCCCTCTGGCGTGGACAATCCGTCAGCACCTGG
 GGAGGAGCCCTGCGCTACCAGCAAGGGAAGATGTCATCCTTGAAGAGGCTCCAGCTCTGCAGATCTGCT
 TCACCAACACCAAGTCCCACGAAGCACCAGGCCCTCGTGGCTGGCGTCAAGAAGCAGGCTAATCAAGTT
 CCCTGAGATCATGGCCCCGCTCCTGACATCAATTGACGCAATCTCCCTGGAGTGTGAGCGCGTGTGGGA
 GAGATGGCGGCCGACCAGTCCCAGAACAGTACCTTGTCTAGAAAGAGCTAATGGACATGAACCAGCACC
 ATCTGAATGCCCTTGGTGTGGGCCACGCTCCCTGGACCAGCTCTGTGAGTAAACAGCAGCACATGGACT
 GCACAGCAAGCTGACTGGCGCAGGCGGGCGGCTGTGGCATCACCTCCTGAAGCCAGGCTAGAGCGA
 GCAAAAGTGGAGGCCCAAGCAGGCCCTGACCGGCTGCGGGTTTGACTGCTGGGAGACCAGCATTGGAG
 CGCCTGGGTCTCCATGCACTCAGCCACCTCCATAGAGGACCCTGTCCGACAAGCCCTGGGCTCCGCT
 GGGTCTCCATGCACTCAGCCACCTCCATAGAGGACCCTGTCCGACAAGCCCTGGGCTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR202119 representing NM_031063
 Red=Cloning site Green=Tags(s)


MLSEVLLVSAPGKVIHGEHAVVHGKVALAVALNLRFTLVLRPQSNQKVSLLNPNVGIKQVVDVATLQLL
 DTGFLEQGDVPAPTLEQLEKLEKLVAGLPRDCVGNESLLAFLYLILAICRKQRTLPPLDIMVWSELPPG
 AGLGSSAAYSVCVAAALLTACEEVTNPLKDRGSIQSWPEEDLKSINKWAYEGERVIHGNPSGVDNSVSTW
 GGALRYQQGKMSSLRPALQIILLTNTKVPSTKALVAGVRSRLIKFPEIMAPLLTSDAISLECERVLG
 EMAAAPVPEQYLVEELMDMNQHHLNALGVGHASLDQLCQVTAHGLHSLKLTGAGGGCGGITLLKPGLER
 AKVEAAKQALTGCGFDCWETSIGAPGVSMHSATSIEDPVQRQALGLRLGSPCTQPPP*RTLSDKPWA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_031063

ORF Size:	1185 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	1723 bp
RefSeq ORF:	1188 bp
Locus ID:	81727
UniProt ID:	P17256
Cytogenetics:	12q16
MW:	42 kDa
Gene Summary:	kinase that catalyzes the ATP-dependent phosphorylation of mevalonic acid to form mevalonate 5-phosphate; involved in the cholesterol biosynthetic pathway [RGD, Feb 2006]