

## Product datasheet for **RR202154**

### Uchl1 (NM\_017237) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Uchl1 (NM\_017237) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Uchl1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR202154 representing NM\_017237  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGCTGAAACCGATGGAGATTAACCCCGAGATGCTGAACAAAGTGTGGCCAAGCTGGGGTGCGCC  
GGCAGTGGCGCTTTGCCGACGTGCTAGGGCTGGAGGAGGAGACTCTGGGCTCAGTGCCATCTCCTGCCTG  
CGCCCTGCTGCTGTTTCCCCTCACGGCCAGCATGAAAACCTCAGGAAAAACAAATGAGGAACTG  
AAGGGACAAGAAGTTAGCCCTAAAGTTTACTTCATGAAGCAGACCATCGGAACTCCTGTGGTACCATTG  
GGCTGATCCACGCACTGGCCAATAACCAAGACAAGCTGGAATTTGAGGATGGATCAGTCCTGAAACAGTT  
TCTGTCTGAAACGGAGAAGTTGTCCCCTGAAGACAGAGCCAAGTGTTCGAGAAGAACGAGGCCATTGAG  
GCAGCCCATGACTCCGTGGCCAGGAGGGCCAGTGCCGGGTAGACGACAAAGTGAATTTCCATTTTATCC  
TGTTCAATAATGTGGACGGCCACCTCTACGAGCTCGATGGGCGAATGCCTTTCCCGGTGAACCATGGCGC  
CAGTTCAGAGGACTCTGCTGCAGGATGCCGCAAGGTCTGCAGAGAATTCAGTGCAGGAGCAGGGA  
GAGGTCCGCTTCCGCGAGTGGCTCTCTGCAAAGCAGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR202154 representing NM\_017237  
Red=Cloning site Green=Tags(s)

MQLKPMIINPEMLNKVLAKLGVAGQWRFADVLGLEEETLGSVSPACALLLLFPLTAQHENFRKKQIEEL  
KGQEVSPKVYFMKQITIGNSCGTIGLIHAVANNQDKLEFEDGSVLKQFLSETEKLSPEDRAKCFEKNEAIQ  
AAHDSVAQEGQCRVDDKVNHFHILFNNVDGHLIELDGRMPFPVNHGASSEDSLLQDAAKVCREFTEREQG  
EVRFSAVALCKAA


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_017237

ORF Size: 669 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_017237.3</a></u> , <u><a href="#">NP_058933.2</a></u>
<b>RefSeq Size:</b>	1058 bp
<b>RefSeq ORF:</b>	672 bp
<b>Locus ID:</b>	29545
<b>UniProt ID:</b>	<u><a href="#">Q00981</a></u>
<b>Cytogenetics:</b>	14p11
<b>MW:</b>	24.8 kDa
<b>Gene Summary:</b>	catalyzes the hydrolysis of ubiquitin carboxy-terminal thiolesters to form ubiquitin and a thiol; may play a role in neuropathic pain [RGD, Feb 2006]