

Product datasheet for **RR202188L3V**

Gde1 (NM_032615) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gde1 (NM_032615) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gde1
Synonyms:	Mir16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032615
ORF Size:	993 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202188).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032615.2 , NP_116004.2
RefSeq Size:	1555 bp
RefSeq ORF:	996 bp
Locus ID:	60418
UniProt ID:	Q9JL55
Cytogenetics:	1q35
Gene Summary:	a putative membrane glycerophosphodiester phosphodiesterase; may mediate lipid metabolism and G protein signaling [RGD, Feb 2006]



[View online »](#)