

## Product datasheet for **RR202190L3V**

### Tram1 (NM\_001007701) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tram1 (NM_001007701) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tram1
Synonyms:	MGC95042
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007701
ORF Size:	1122 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202190).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001007701.1</a> , <a href="#">NP_001007702.1</a>
RefSeq Size:	2778 bp
RefSeq ORF:	1125 bp
Locus ID:	312903
UniProt ID:	<a href="#">Q5XI41</a>
Cytogenetics:	5q11
Gene Summary:	Stimulatory or required for the translocation of secretory proteins across the ER membrane. [UniProtKB/Swiss-Prot Function]



[View online »](#)