

Product datasheet for **RR202219L3V**

Osmr (NM_001005384) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Osmr (NM_001005384) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Osmr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001005384
ORF Size:	2886 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202219).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001005384.1
RefSeq Size:	2889 bp
RefSeq ORF:	2889 bp
Locus ID:	310132
UniProt ID:	Q65Z14
Cytogenetics:	2q16
Gene Summary:	mouse homolog binds cytokine signalling subunits c12 protein and gp130 and acts as a receptor for the cytokine oncostatin M [RGD, Feb 2006]



[View online »](#)