

Product datasheet for RR202223

P2rx6 (NM_012721) Rat Tagged ORF Clone

Product data:

| | |
|--------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | P2rx6 (NM_012721) Rat Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | P2rx6 |
| Synonyms: | P2rx1; P2x6; Rnap2x6 |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| ORF Nucleotide Sequence: | >RR202223 representing NM_012721 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGCAGTTGCAGCCAGCGGAACTGGGAGCATGGCTTCTGCAGTCGCAGCAGCACTAGTGAGCTGGGGT
TTCTGGATTACAAGACAGAGAAATATGTGATGACCAGAACTGCTGGGTGGGCATCTCCAGAGGCTGCT
GCAGCTCGGAGTGGTAGTCTACGTGATTGGGTGGGCCCTCTGGCCAAAAAAGGCTACCAGGAATGGGAC
ATGGACCCCGAGATTTCTGTCATCACAAACTCAAGGGGTTTCTGTAACCCAGGTTAAGGAACTGGAGA
AACGGCTGTGGGATGTGGCTGACTTCGTGAGGCCATCACAGGGAGAGAATGTGTTCTTCTGGTAACCA
CTTCTTGAGACACCAGCTCAAGTCCAGGGCAGATGCCAGAGCATCCTTCTGTTCTCTGGCTAACTGC
TGGGCTGATGAGGACTGCTCTGAAGGGGAGATGGGGACATACAGCCACGGCATAAAGACTGGTCAGTGTG
TGGCGTTCAACGGAACCCACAGGACCTGTGAGATCTGGAGCTGGTGGCCAGTGGAGAGTAGTGTGTGCC
CAGGAAACCCCTGCTAGCTCAGGCCAAGAATTACGCTCTTCATCAAAAATACTGTAACCTTCAACAAG
TTCAACTTCTCCAGAACCAATGCCTTAGATACCTGGGACAACACCTATTTCAAGTACTGTCTCTACGATT
CACTCTCCAGTCCGTACTGCCAGTGTTCGCATCGGGGACCTTGTGGCTAGACTGGAGGGGACTTTGA
GGACCTTGCATTGCTGGGTGGTGTGTGGGCATCAACATCCACTGGGATTGCAACCTGGACACGAAAGGC
TCTGACTGCTCTCCCAGTACTCCTCCAGCTGCAGGAGAGGGGGTACAATTTCCAGGACAGCTAATTACT
GGTGGGACGCTCAGGCTGGAGTCCCAGGACCTGCTCAAACCTATGGGATCCGCTTTGACATTCTGT
CACTGGGACAGGAAAGTTTGCCTACTCATCCCTACGGCCATCACAGTGGGCACAGGAGCAGCTTGCTG
GGCATGGTCACTTCTCTGTGATCTGTGCTACTGTACGTGGATAGAGAGGCCGTTTCTACTGGAGGA
CCAAGTATGAAGAGGCAAGAGCCCCAAAGGCAACTACCAACAGTGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR202223 representing NM_012721
 Red=Cloning site Green=Tags(s)


MQLQPAGTGMASAVAAALVSWGFLDYKTEKYVMTRNCWVGISQRLQLGVVVYVIGWALLAKKGYQEWDMDPQISVITKLKGVSVTQVKELEKRLWDVADFVRPSQGENVFFLVTFNFLTVAQVQGRCPPEHPSVPLANCWADEDCEGEMGTYSHGIKTGQCVAFNGTHRTCEIWSWCPVESSAVPRKPLLAQAKNFTLF IKNTVTFNKFNF SRTNALDTWDNTYFKYCLYDSLSSPYCPVFRIGDLVAMTGGDFEDLALLGGAVGINIHWCNLDTKGSDCSPQYSFQLQERGFNFRTANYWAAASGVESRSLKLKYGIRFDILVTGQAGKFALIPTAITVGTGAAWLGMVTFLCDLLLLYVDREAGFYWRTKYEEARAPKATTNSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_012721

| | |
|-------------------------------|---|
| ORF Size: | 1167 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_012721.2 , NP_036853.2 |
| RefSeq Size: | 2331 bp |
| RefSeq ORF: | 1170 bp |
| Locus ID: | 25041 |
| UniProt ID: | P51579 |
| Cytogenetics: | 11q23 |
| MW: | 43.4 kDa |
| Gene Summary: | acts as an ATP-gated ion channel [RGD, Feb 2006] |