

Product datasheet for **RR202238L3V**

Apmmap (NM_001034003) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apmmap (NM_001034003) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Apmmap
Synonyms:	Ab2-305; RGD1308874
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001034003
ORF Size:	1128 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202238).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001034003.1 , NP_001029175.1
RefSeq Size:	2099 bp
RefSeq ORF:	1131 bp
Locus ID:	366227
UniProt ID:	Q7TP48
Cytogenetics:	3q41
Gene Summary:	Exhibits strong arylesterase activity with beta-naphthyl acetate and phenyl acetate. May play a role in adipocyte differentiation (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)