

Product datasheet for **RR202249L4V**

Bysl (NM_182674) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bysl (NM_182674) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bysl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182674
ORF Size:	1308 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182674.3 , NP_872615.3
RefSeq Size:	1696 bp
RefSeq ORF:	1311 bp
Locus ID:	359727
UniProt ID:	Q80WL2
Cytogenetics:	9q12
Gene Summary:	Required for processing of 20S pre-rRNA precursor and biogenesis of 40S ribosomal subunits.[UniProtKB/Swiss-Prot Function]



[View online »](#)