

## Product datasheet for **RR202284L4V**

### Clpb (NM\_022947) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Clpb (NM_022947) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Clpb
Synonyms:	Skd3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022947
ORF Size:	2031 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022947.1</a> , <a href="#">NP_075236.1</a>
RefSeq Size:	2138 bp
RefSeq ORF:	2034 bp
Locus ID:	65041
UniProt ID:	<a href="#">Q9WTT2</a>
Cytogenetics:	1q32
Gene Summary:	mouse homolog is a putative ATPase and member of the Clp/HSP104 family [RGD, Feb 2006]



[View online »](#)