

Product datasheet for **RR202304**

Ikzf1 (NM_001107237) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Ikzf1 (NM_001107237) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Ikzf1
 Synonyms: RGD1562979
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR202304 representing NM_001107237
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGATGTCGATGAGGGTCAAGACATGTCCAAGTTTCAGGAAAGGAGAGCCCTCCAGTCAGTGACACTC
 CGGATGAAGGGGATGAACCCATGCCTGTCCCGGAGGACCTGTCCACCTCTGGAGCACAGCAGAATCCAA
 GAGTGATCGAGGCATGGGTGAACGGCCTTTCCAGTGCAACCAGTGTGGGGCTCCTTTACCCAGAAAGGC
 AACCTACTGCGGCACATCAAGCTGCACTCGGGCGAGAAGCCCTTCAAATGCCATCTCTGCAACTATGCCT
 GCCCGCGAGGGACGCCCTCACCGGCCACCTAAGGACGCACTCCGTTGGTAAGCCTCACAAATGTGGATA
 TTGTGGCCGGAGCTATAAACAGCGAAGCTTTTAGAAGAACATAAAGAGCGATGCCACAACACTACTGGAA
 AGCATGGGCCCTCCCGGCATGTACCCAGTCATTAAGGAAGAACTAACCACAGCGAAATGGCAGAGGACC
 TGTGCAAGATAGGAGCAGAGAGGTCCCTTGTCTGGACAGGCTGGCCAGCAACGTCGCCAAACGTAAGAG
 CTCTATGCCTCAGAAATTTCTTGAGACAAGTGCCTGTGACAGATGCCCTATGACAGTGCCAACATGAG
 AAGGATGAGATGATGACATCCACGTGATGGACCAGGCCATCAACAATGCCATTAACACTCTGGGGCTG
 AGTCCCTGCGCCGTTGGTGCAGACACCCCTGGTAGCTCCGAGGTGGTCCCGGCATCAGTTCCCATGTA
 CCAGCTGCACAAGCCCCCTCAGACGGGCCACCGGTCCAACCACTCAGCACAGGACAGCGCCGTGGAG
 AACTTGCTGCTGTCTCAAGGCCAAGTCCGTGTATCAGAGCGGGAGGCCTCCCCAAGCAACAGCTGCC
 AAGACTCCACGGATACAGAGAGCAACCGGGAGGAACAGCGCAGCGGTCTCATCTACCTAACCAACCAT
 CACTCCGCATGCGCGCAACGGGCTGGCTCTCAAGGAGGAGCAGCGCGCCTACGAGGTGCTGAGGGCAGCC
 TCAGAGAACTCGCAGGACGCCTCCGTGTGGTCAAGGAGGAGTGGTGAAGGTGTACAAGTGGC
 AACACTGCCGCTGCTTCTCTGGATCATGTATACACCATCCACATGGGCTGCCATGGCTTCCGGGA
 TCCCTTTGAGTGAACATGTGTGGTTATCACAGCCAGGACAGGTACGAGTTCTCTCCACATCACAGG
 GGGGAGCATCGTACCACCTGAGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR202304 representing NM_001107237
 Red=Cloning site Green=Tags(s)



MDVDEGQDMSQVSGKESPPVSDTPDEGDEPMPVPEDLSTSGAQQNSKSDRGMGERPFQCNQCGASFQKG
 NLLRHIKLSHGEKPFKCHLCNYACRRRDALTGHLRTHSVGKPHKCGYCGRSYKQRSSLEEHKERCHNYLE
 SMGLPGMYPVIKEETNHSEMAEDLCKIGAERSLVLDRLASNAVKRKSMPQKFLGDKCLSDMPYDSANYE
 KDEMMTSHVMDQAINNAINYLGAESLRPLVQTPPGSSEVVPVISSMYQLHKPPSDGPPRSNHSQAQDSAVE
 NLLLLSKAKSVSSEREASPSNSCQDSTDTESNAEEQRSGLIYLTNHI TPHARNGLALKEEQRAYEVLRAA
 SENSQDAFRVVSTSGEQLKVKCEHCRVLF LDHVMYTIHMGGCHGRDPFECNMGYHSQDRYEFSSHITR
 GEHRYHLS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001107237

ORF Size:	1284 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001107237.1 , NP_001100707.1
RefSeq Size:	2023 bp
RefSeq ORF:	1287 bp
Locus ID:	305501
Cytogenetics:	14q21
MW:	47.9 kDa