

Product datasheet for **RR202388L3V**

Nanp (NM_001009409) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nanp (NM_001009409) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nanp
Synonyms:	Hdhd4; RGD1306009
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001009409
ORF Size:	744 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202388).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009409.1 , NP_001009409.1
RefSeq Size:	1335 bp
RefSeq ORF:	747 bp
Locus ID:	311530
UniProt ID:	Q5M969
Cytogenetics:	3q41
Gene Summary:	The protein encoded by this gene is an enzyme involved in the synthesis of N-acetylneuraminate (Neu5Ac), the main form of sialic acid. [provided by RefSeq, Jul 2008]



[View online »](#)