

## Product datasheet for **RR202397L4V**

### **Pcdha5 (NM\_199505) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pcdha5 (NM_199505) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pcdha5
Synonyms:	rCNRv05
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199505
ORF Size:	2802 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR202397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199505.4</a> , <a href="#">NP_955799.1</a>
RefSeq Size:	5385 bp
RefSeq ORF:	2805 bp
Locus ID:	393087
Cytogenetics:	18p11
Gene Summary:	member of a family of synaptic cadherins and Reelin receptors, integral plasma membrane proteins that may play a critical role in maintaining specific cell-cell interactions in the brain [RGD, Feb 2006]



[View online »](#)